

## StarLab Release Notes

---

### Version 3.93 (.02) Released 23-December-2024

#### *New features:*

- The COM object now supports wireless connections to Ariel

#### *Updates:*

- COM objects version 10.8
- Updated "Ophir User Commands"

#### *Fixes:*

- Allow DOSAGE mode for StarLite, with PD300-RM

#### *Known Issues:*

- Time-out when searching for non-existing device on Bluetooth or Ethernet can take up to a minute!

## Version 3.92 (.03) Released 11-July-2024

### *New features:*

- Added PD300RM support for Centauri on StarLab
- Enhancement: Add auto scale option for line graphs
- Support connecting to Centauri via Ethernet
- Extended "Time Frame" of graphs to be up to 99:59 hours

### *Updates:*

- Update "Ophir User Commands"
- COM objects version 10.7

### *Fixes:*

- StarLab previous version didn't run on Win 7
- When attempting to open a second instance of StarLab there was an assertion failure

### *Known Issues:*

- Time-out when searching for non-existing device on Bluetooth or Ethernet can take up to a minute!